

Feather LIFT™



What is “Feather-Lift” ?

The “Feather-Lift” is a non-surgical approach to facial lifting utilizing the patented “Aptos” threads.

What are “Aptos” threads ?

“Aptos” threads are made of a monofilament material called “polypropylene”. These “Aptos” threads are designed with directional “cogs”. Once inserted under the skin, the “Aptos” thread cogs will form a support structure for the tissues of the face, while preventing movement.



Who is the best candidate ?

The best candidate is a patient with:

- ⇒ Ptosis of the soft tissues of the face and neck.
- ⇒ Premature aging and sun damaged face.
- ⇒ Weakly pronounced aesthetic contours.
- ⇒ Flaccid, flat face


Is it hard to perform ?

The procedure is quick and easy to perform. It is done under local anesthesia and if desired, light IV sedation, the patient remaining awake and alert. The “Aptos” threads are inserted using a needle in the subcutaneous fat, at the required depth, along the contours that were previously marked.



How does the “Feather-Lift” differ from conventional face lift ?

- ⇒ Uniform gathering and even lifting of the soft tissues for a new volumetric contour of the face.
- ⇒ Fast recovery, patient may return to work in 1 week.
- ⇒ No scars on the skin.



Suture with directional cogs.